

Mini Stick max

Coaxial Microintroducer Kit

STARTING HAS NEVER BEEN EASIER.

AngioDynamics is committed to providing physicians the best tools possible for the best patient outcomes, no matter how big or small the device. Designed to start each procedure with ease and efficiency, the AngioDynamics Mini Stick MAX Coaxial Microintroducer Kits incorporate proprietary designs and features that aid in simplifying set up and gaining vascular access.

NEEDLE OPTIONS

4 cm and 7 cm length
available with or without
echogenic tip

SHEATH DILATOR OPTIONS

10 cm and 15 cm
sheath dilator options

Containment Clip

designed to keep all components
together in a sterile field



THREE DIFFERENT 0.018" WIRE CONFIGURATIONS

- Stainless/Tungsten
- Nitinol/Tungsten
- Stainless Steel/Stainless Steel

Translucent hub
enables visualization
of flash back

Interlocking coaxial Design

in 4F and 5F with stiff
or standard dilator




angiodynamics

GL/VA/BR/1973 Rev 01 08/2023

Tapered dilator
for a smooth transition at the tip

Proprietary coaxial design
for enhanced ease of insertion



Dimpling on distal tip
for enhanced echogenicity

Black bevel indicator
provides easy indication of
bevel orientation

KIT COMPONENTS							MINI STICK MAX COAXIAL MICROINTRODUCER KITS ^{1,2}				BOX QTY
0.018" WIRE		21G NEEDLE				SHEATH	STANDARD		STIFFENED		
CORE/TIP MATERIAL	LENGTH (CM)	ECHOGENIC TIP	LENGTH (CM)	BEVEL	HUB	LENGTH (CM)	4F	5F	4F	5F	
NITINOL / TUNGSTEN	45	YES	4	B	GREEN	10	H96545771	H96545775	H96545781	H96545785	1
			7	B	GREEN	10	H96545750	H96545756	H96545753	H96545759	1
		NO	4	B	GREEN	10	H96545769	H96545773	H96545778	H96545783	1
				B	GREEN	15				H96545786	1
			7	B	GREEN	10	H96545763	H96545767	H96545762	H96545766	1
				B	GREEN	15	-	-	-	H96545776	1
STAINLESS STEEL / STAINLESS STEEL	45	YES	4	B	GREEN	10	H96545770	H96545774	H96545779	H96545784	1
			7	B	GREEN	10	H96545748	H96545754	H96545751	H96545757	1
		NO	4	B	GREEN	10	H96545768	H96545772	H96545777	H96545782	1
			7	B	GREEN	10	H96545761	H96545765	H96545760	H96545764	1
STAINLESS STEEL / TUNGSTEN	45	YES	7	B	GREEN	10	H96545749	H96545755	H96545752	H96545758	1

Each kit includes: one .018" guidewire, one 10 cm sheath dilator/introducer and one 21 gauge needle† 15 cm sheath

FOOTNOTES:

¹Not all products may be available for sale in your market. Please consult with your local AngioDynamics Partner for product availability and additional product information.

²Low volume configurations may be built to order and require a greater lead time.

IMPORTANT RISK INFORMATION

MINI STICK MAX COAXIAL MICROINTRODUCER KIT

INDICATION FOR USE: The coaxial microintroducer kit is used for the percutaneous introduction of a guidewire into the vascular system.

CONTRAINDICATIONS: None known.

Refer to Directions for Use and/or User Manual provided with the product for complete Instructions, Warnings, Precautions, Possible Adverse Effects and Contraindications prior to use of the product.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

USA > 14 Plaza Drive, Latham, NY 12110 > Tel: 800-772-6446 or 518-798-1215 > Fax: 518-798-1360
International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
Tel: +31 (0)20 753 2949 > Fax: +31 (0)20 753 2939

www.angiodynamics.com



*AngioDynamics, the AngioDynamics logo and Mini Stick are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. All other trademarks are property of their respective owners. © 2023 AngioDynamics, Inc. GL/VA/BR/1973 Rev 01 08/2023